


Product Datasheet

JUND antibody (orb176721)

Description	Goat polyclonal antibody to JUND.
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	JUND
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 40kDa band observed in nuclear lysates of cell line NIH3T3 (calculated MW of 35.2kDa according to NP_005345.3). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Short and Pfarr, J Biol Chem. 2002
Protein Sequence	CQLLPQHQPAY
Clonality	Polyclonal
MW	35.2
Entrez	3727
Dilution Range	ELISA: 1:64000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH3T3 nuclear ...