


Product Datasheet

Axotrophin antibody (orb176725)

Description	Goat polyclonal antibody to Axotrophin.
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MARCH7 / Axotrophin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 75kDa band observed in Mouse Skeletal Muscle lysates (calculated MW of 76.6kDa according to NP_065600.1. Recommended concentration: 1-3µg/ml.
Protein Sequence	RTLQAHMEDLETSED
Clonality	Polyclonal
MW	78.1
Entrez	64844
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Skeletal ...