


Product Datasheet

MEX3C antibody (orb176726)

Description	Goat polyclonal antibody to MEX3C.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	MEX3C (aa541-554)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 60+36kDa bands observed in Mouse Testis lysates, and these bands reduced after knock-down by gene trap (calculated MW of 69.4kDa according to NP_057710.3). Recommended concentration: 0.2-0.5µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Kidney. Recommended concentration: 5µg/ml.
Protein Sequence	TPRLSPTFPESIEH
Clonality	Polyclonal
MW	69.4
Entrez	51320
Dilution Range	ELISA: 1:32000, WB: 0.2-0.5 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Testis US...