


Product Datasheet

OAS1 antibody (orb176728)

Description	Goat polyclonal antibody to OAS1.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	OAS1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50kDa band observed in lysates of cell line HepG2 (calculated MW of 46.0kDa according to NP_058132.2). Recommended concentration: 1-3µg/ml.
Protein Sequence	THEYPHFSHRPST
Clonality	Polyclonal
MW	46.0; 41.7; 47.4
Entrez	4938
Dilution Range	P-ELISA: 1:128000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 lysate US...