


Product Datasheet

Syndecan 1 antibody (orb176733)

Description	Goat polyclonal antibody to Syndecan 1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	syndecan 1 / CD138
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 37kDa band observed in preliminary testing of HepG2 cell lysate after 1ug/ml of antibody staining (calculated MW of 32.5kDa according to NP_002988.3). Primary incubation 1 hour at room temperature.
Protein Sequence	REQEATPRPRETTQ
Clonality	Polyclonal
MW	32.5
Entrez	6382
Dilution Range	ELISA: 1:2000, WB: 1ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 lysate us...