


Product Datasheet

syntaxin 11 antibody (orb176734)

Description	Goat polyclonal antibody to syntaxin 11.
Species/Host	Goat
Reactivity	Canine, Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	syntaxin 11
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 36kDa band observed in Human Placenta lysates (calculated MW of 33.2kDa according to NP_003755.2). This molecular weight is routinely observed by other sources. Recommended concentration: 1-3µg/ml.
Protein Sequence	KQADTLNVIELNVQK
Clonality	Polyclonal
MW	33.2
Entrez	8676
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Placenta lysate...