


Product Datasheet

TTLL4 antibody (orb176737)

Description	Goat polyclonal antibody to TTLL4.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TTLL4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 150kDa band observed in lysates of cell line MDCK (calculated MW of 133kDa according to NP_055455.3). Recommended concentration: 1-3µg/ml.
Protein Sequence	RSSYRQKPYQQLE
Clonality	Polyclonal
MW	133
Entrez	9654
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of MDCK lysate using TTLL4 antibody.