

Product Datasheet

Anti Mouse Ep-CAM (CD326 / TACSTD1) Mab [323/A3] [Out of stock] (orb1768089)

Description	Anti Mouse Ep-CAM (CD326 / TACSTD1) Mab [323/A3] [Out of stock]
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Note	For research use only
Isotype	Mouse
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.