

Product Datasheet

Anti Human Stat3 (Phospho-Y705) Pab Rabbit [Out of stock] (orb1769932)

Description	Anti Human Stat3 (Phospho-Y705) Pab Rabbit [Out of stock]
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Note	For research use only
Isotype	Rabbit
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.