

Product Datasheet

Anti Human MEK-3 (Phospho-S218) Pab Rabbit [Out of stock] (orb1772782)

Description	Anti Human MEK-3 (Phospho-S218) Pab Rabbit [Out of stock]
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Note	For research use only
Isotype	Rabbit
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.