

Product Datasheet

Anti Human EGFR (Phospho-Y1069) Pab Rabbit [Out of stock] (orb1776791)

Description	Anti Human EGFR (Phospho-Y1069) Pab Rabbit [Out of stock]
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Note	For research use only
Isotype	Rabbit
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.