

Product Datasheet

Anti Human CD155 (TAGE4 / PVR) Mab [SK II 4.3] [Out of stock] (orb1779102)

Description	Anti Human CD155 (TAGE4 / PVR) Mab [SK II 4.3] [Out of stock]
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, IP, WB
Note	For research use only
Isotype	Mouse
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.