

Product Datasheet

Anti Human Apolipoprotein B (ApoB) Mab [Out of stock] (orb1780425)

Description	Anti Human Apolipoprotein B (ApoB) Mab [Out of stock]
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Human ApoB native protein
Preservatives	Prepared as a sterile filtered solution in phosphate buffer, pH 7.2
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	For short term storage, store at 4oC for up to one month. For long term storage, store at -20oC or -80oC. May be made into aliquots. Avoid repeated freeze thaw cycles
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Affinity purified antibody, > 99% pure by SDS PAGE
Expiration Date	12 months from date of receipt.