

Product Datasheet

Anti FOXP3 (Forkhead Box Protein P3) / Scurfin Mab [FXP3/197] [Out of stock] (orb1781825)

Description	Anti FOXP3 (Forkhead Box Protein P3) / Scurfin Mab [FXP3/197] [Out of stock]
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Note	For research use only
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.