



Product Datasheet

MYO1F antibody (orb178309)

Description	Rabbit polyclonal antibody to MYO1F
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 384-592aa of human MYO1F protein
Target	MYO1F
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Concentration	0.5mg/ml
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	124kDa
Dilution Range	WB: 1:800-1:2,000, ELISA: 1:20,000-1:80,000, IHC-P: 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lys...



Immunohistochemical staining of fetal st...