


Product Datasheet

Human VEGF 165 protein (orb178366)

Description	Human VEGF 165 protein
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human VEGF 165
Form/Appearance	0.14mol/L NaCl, 0.027mol/L KCl, 0.01mol/L Na ₂ PO ₄ , 0.018mol/L KH ₂ PO ₄
Concentration	- 100 µg (in 100 µl): 1 mg/ml- 500 µg (in 0.5 µl): 1 mg/ml
Storage	Store at 4°C for short term and -20°C for long term storage
Note	For research use only
Application notes	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Purity	> 85% as determined by SDS-PAGE
Uniprot ID	P15692
NCBI	003376, 003367
Entrez	7422
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of human vascular endo...