



Product Datasheet

GOT2 antibody (orb178369)

Description	Goat polyclonal antibody to GOT2.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GOT2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 45kDa band observed in lysates of cell line HepG2, and of Mouse Brain and Heart (calculated MW of 47.5kDa according to human NP_002071.2 and 47.4kDa according to mouse NP_034455.1). Recommended concentration: 0.3-1µg/ml.
Protein Sequence	QGYRYDPKTCGFD
Clonality	Polyclonal
MW	47.5
Entrez	2806
Dilution Range	ELISA: 1:8000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 lysate US...



Western blot analysis of Mouse Brain (La...