

Product Datasheet

alpha 2 Macroglobulin antibody (orb178373)

Description Goat polyclonal antibody to alpha 2 Macroglobulin

Species/Host Goat

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Target alpha-2-macroglobulin

Preservatives Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage Aliquot and store at -20°C. Minimize freezing and thawing.

Note For research use only

Application notes WB: Approx 200kDa band observed in Human Albumin-depleted Serum lysates (calculated MW of 163kDa according to NP_000005.2). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 0.5-2µg/ml. Primary incubation was 1 hou. IHC: Paraffin embedded Human Placenta and Liver. Recommended concentration: 5µg/ml.

Protein Sequence TSLTVRVNYKDRSP


Clonality Polyclonal

MW 163

Entrez [2](#)

Dilution Range ELISA: 1:32000, WB: 0.5-2 µg/ml, IHC-P: 5 µg/ml

Expiration Date 12 months from date of receipt.



Western blot analysis of Human Albumin-d...