


Product Datasheet

OAS2 antibody (orb178374)

Description	Goat polyclonal antibody to OAS2.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	OAS2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 85 +19kDa bands observed in lysates of cell line Jurkat (calculated MW of 82.4kDa according to NP_058197.2 and of 19.5kDa according to NP_001027903.1). Recommended concentration: 1-3µg/ml.
Protein Sequence	RKTVLRGNSDGT
Clonality	Polyclonal
MW	82.4; 78.8; 19.5
Entrez	4939
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat lysate u...