



Product Datasheet

OAS2 antibody (orb178375)

Description	Goat polyclonal antibody to OAS2.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	OAS2 (aa357-371)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 70kDa band observed in lysates of cell lines Daudi, Jurkat and K562 (calculated MW of 78.8kDa according to NP_002526.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Lin et al, J Immunol. 20. IHC: Paraffin embedded Human Breast. Recommended concentration: 5µg/ml.
Protein Sequence	QPNKCFLEQIDSAVN
Clonality	Polyclonal
MW	82.4; 78.8; 19.5
Entrez	4939
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Daudi (Lane 1),...



Immunohistochemical staining of Human Br...