



Product Datasheet

LDHB antibody (orb178383)

Description	Goat polyclonal antibody to LDHB.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	LDHB (aa15-25)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 35kDa band observed in Human Heart and Kidney lysates (calculated MW of 36.6kDa according to NP_002291.1). Recommended concentration: 0.3-1µg/ml. IHC: Paraffin embedded Human Kidney. Recommended concentration: 5µg/ml.
Protein Sequence	EEATVPNNKIT
Clonality	Polyclonal
MW	36.6
Entrez	3945
Dilution Range	ELISA: 1:4000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart (La...)



Immunohistochemical staining of Human Ki...