


Product Datasheet

SERPINB1 antibody (orb178384)

Description	Goat Polyclonal antibody to SERPINB1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	SERPINB1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 43kDa band observed in lysates of cell line A549 (calculated MW of 42.7kDa according to NP_109591.1). Recommended concentration: 1-3µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Spleen (neutrophils). Recommended concentration: 5µg/ml.
Protein Sequence	EVHSRFQSLNADINKR
Clonality	Polyclonal
MW	42.7
Entrez	1992
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 lysate using...