

Product Datasheet

CD48 Antibody (BSA Free) (orb1784564)

Description	CD48 Antibody (BSA Free)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	residues 23-216 (ECD)
Target	CD48
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, rAlbumin-Free
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P18181
Expiration Date	12 months from date of receipt.