



Product Datasheet

GFP Antibody (orb1784590)


Description	GFP Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Green Fluorescent Protein
Target	GFP
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	50 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Dilution Range	WB: 1:1000, IP: 20 µl/mg lysate, ICC: 1:100 to 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE cell sections
Expiration Date	12 months from date of receipt.



Detection of GFP by western blot. Samples...



Detection of GFP by western blot. Samples...



Detection of GFP by western blot of immu...