


## Product Datasheet

### LTA4(Leukotriene A4)ELISA Kit (orb1784592)

|                        |                               |
|------------------------|-------------------------------|
| <b>Description</b>     | LTA4(Leukotriene A4)ELISA Kit |
| <b>Reactivity</b>      | All                           |
| <b>Range</b>           | 15.625-1000pg/ml              |
| <b>Target</b>          | LTA4                          |
| <b>Storage</b>         | 2-8 °C for 6 months           |
| <b>Note</b>            | For research use only         |
| <b>Sensitivity</b>     | 9.375ng/ml                    |
| <b>Expiration Date</b> | Please enquire.               |



LTA4(Leukotriene A4) ELISA Kit