


Product Datasheet

Rat PFKFB3(6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 3) ELISA Kit (orb1784609)

Description	Rat PFKFB3(6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 3) ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	PFKFB3
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat PFKFB3(6-phosphofructo-2-kinase/fruc...