


Product Datasheet

Rat Anti SARS-CoV-2 (COVID) Spike RBD Alpha B.1.1.7 Variant(UK Variant) IgG ELISA Kit (orb1784614)

Description	Rat Anti SARS-CoV-2 (COVID) Spike RBD Alpha B.1.1.7 Variant(UK Variant) IgG ELISA Kit
Reactivity	Rat
Range	1.56-100ng/ml
Target	Anti SARS-CoV-2 UK Variant IgG
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Rat Anti SARS-CoV-2 (COVID) Spike RBD Al...