


Product Datasheet

Rat MAS-R(Mas-related G protein-coupled receptor b)ELISA Kit (orb1784620)

Description	Rat MAS-R(Mas-related G protein-coupled receptor b)ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	MAS-R
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	W8W3K0
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat MAS-R(Mas-related G protein-coupled
...
...