


Product Datasheet

Monkey RF-IgG(Rheumatoid Factor IgG) ELISA Kit (orb1784645)

Description	Monkey RF-IgG(Rheumatoid Factor IgG) ELISA Kit
Reactivity	Monkey
Range	3.125-200IU/ml
Target	RF-IgG
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	1.875IU/ml
Expiration Date	Please enquire.



Monkey RF-IgG(Rheumatoid Factor IgG) ELISA Kit