


Product Datasheet

Mouse Calcrl (Calcitonin gene-related peptide type 1 receptor) ELISA Kit (orb1784652)

Description	Mouse Calcrl (Calcitonin gene-related peptide type 1 receptor) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Calcrl
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
Calcrl
(Calcitonin
gene-
related
pe...