


Product Datasheet

Mouse Glrx5(Glutaredoxin-related protein 5, mitochondrial) ELISA Kit (orb1784654)

Description	Mouse Glrx5(Glutaredoxin-related protein 5, mitochondrial) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	GlrX5
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
GlrX5(Glutaredoxin-
related protein...