


Product Datasheet

Mouse CCDC25(Coiled-coil domain-containing protein 25)ELISA Kit (orb1784663)

Description	Mouse CCDC25(Coiled-coil domain-containing protein 25)ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	CCDC25
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
CCDC25(Coiled-
coil domain-
containing...