


Product Datasheet

Mouse ISG15(Ubiquitin-like protein ISG15) ELISA Kit (orb1784668)

Description	Mouse ISG15(Ubiquitin-like protein ISG15) ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	ISG15
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse
ISG15(Ubiquitin-
like protein
ISG15...