


Product Datasheet

Human MOTS-c(Mitochondrial-derived peptide MOTS-c)Elisa kit (orb1784671)

Description	Human MOTS-c(Mitochondrial-derived peptide MOTS-c)Elisa kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	MOTS-c
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human MOTS-c(Mitochondrial-derived pepti...