


Product Datasheet

Human anti-ASIC1(Acid-sensing ion channel 1)antibody Elisa kit (orb1784683)

Description	Human anti-ASIC1(Acid-sensing ion channel 1)antibody Elisa kit
Reactivity	Human
Range	3.125-200ng/ml
Target	anti-ASIC1
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	1.875ng/ml
Expiration Date	Please enquire.



Human anti-ASIC1(Acid-sensing ion channel...)