


Product Datasheet

Human Anti SARS-CoV-2 (COVID) Spike RBD Alpha B.1.1.7 Variant(UK Variant) IgA ELISA Kit (orb1784693)

Description	Human Anti SARS-CoV-2 (COVID) Spike RBD Alpha B.1.1.7 Variant(UK Variant) IgA ELISA Kit
Reactivity	Human
Range	0.781-50ng/ml
Target	Anti SARS-CoV-2 UK Variant IgA
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.469ng/ml
Expiration Date	Please enquire.



Human
Anti
SARS-
CoV-2
(COVID)
Spike
RBD ...