

Product Datasheet

Human Anti SARS-CoV-2 (COVID) Spike RBD Alpha B.1.1.7 Variant(UK Variant) IgG ELISA Kit (orb1784695)

Description	Human Anti SARS-CoV-2 (COVID) Spike RBD Alpha B.1.1.7 Variant(UK Variant) IgG ELISA Kit
Reactivity	Human
Range	1.563-100ng/ml
Target	Anti SARS-CoV-2 UK Variant IgG
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.