


Product Datasheet

Human YWHAZ(14-3-3 protein zeta/delta)ELISA Kit (orb1784696)

Description	Human YWHAZ(14-3-3 protein zeta/delta)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	YWHAZ
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
YWHAZ(14-
3-3 protein
zeta/delta)
E...