


Product Datasheet

Human MAS-R(Mas-related G-protein coupled receptor MRG)ELISA Kit (orb1784698)

Description	Human MAS-R(Mas-related G-protein coupled receptor MRG)ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	MAS-R
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P35410
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
MAS-
R(Mas-
related G-
protein
couple...