


Product Datasheet

Human BDH2 (3-hydroxybutyrate dehydrogenase type 2) ELISA Kit (orb1784703)

Description	Human BDH2 (3-hydroxybutyrate dehydrogenase type 2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	BDH2
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human BDH2
(3-
hydroxybutyrate
dehydrogen...