


Product Datasheet

Chicken MASP2(Mannan-binding lectin serine protease 2) ELISA Kit (orb1784709)

Description	Chicken MASP2(Mannan-binding lectin serine protease 2) ELISA Kit
Reactivity	Gallus
Range	31.2-2000pg/ml
Target	MASP2
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Chicken
MASP2(Mannan-
binding lectin
seri...