


Product Datasheet

Canine GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit (orb1784717)

Description	Canine GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit
Reactivity	Canine
Range	15.625-1000pg/ml
Target	GM-CSF
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Canine GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit