

Product Datasheet

anti- TOMM20 antibody (orb1784836)

Description	anti- TOMM20 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	translocase of outer mitochondrial membrane 20 homolog
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2b
Purity	≥95% as determined by SDS-PAGE
MW	16 kDa
Uniprot ID	Q15388
Dilution Range	WB: 1:1000 - 1:10000; IHC: 1:200 - 1:1000; IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.