

Product Datasheet

SARS-CoV-2 NUCLEOPROTEIN (orb1784848)

Description	SARS-CoV-2 NUCLEOPROTEIN
Reactivity	Virus
Tested Applications	ELISA
Preservatives	None present. 20mM Sodium Phosphate25mM Potassium Carbonate pH 10.0150mM Sodium Chloride
Form/Appearance	Purified recombinant protein - liquid
Concentration	Current, batch-specific concentration 1.88 mg/ml
Storage	-70°C (ship frozen)
Note	For research use only
Expiration Date	6 months from date of receipt.