

Product Datasheet

MOUSE ANTI HUMAN GRB10 (orb1784850)

Description	MOUSE ANTI HUMAN GRB10
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	E. coli-derived recombinant protein, aa 1-282 of human GRB10
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	IgG concentration 1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.