

Product Datasheet

RAT ANTI HUMAN RNA POL II CTD (orb1784867)

Description	RAT ANTI HUMAN RNA POL II CTD
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, IF, WB
Immunogen	A synthetic peptide containing the CTD of human RNA polymerase II
Preservatives	0.035% Sodium Azide (NaN ₃)30% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	IgG concentration 1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.