

## Product Datasheet

### RABBIT ANTI HUMAN CTCF (orb1784958)

|                            |                                                                                |
|----------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>         | RABBIT ANTI HUMAN CTCF                                                         |
| <b>Species/Host</b>        | Rabbit                                                                         |
| <b>Reactivity</b>          | Human                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                   |
| <b>Tested Applications</b> | ChIP, IF, IHC-P, WB                                                            |
| <b>Immunogen</b>           | A peptide within the N-terminal region of human CTCF                           |
| <b>Preservatives</b>       | 0.035% Sodium Azide (NaN <sub>3</sub> )30% Glycerol. Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                                          |
| <b>Storage</b>             | +4°C, -20°C if preferred                                                       |
| <b>Note</b>                | For research use only                                                          |
| <b>Isotype</b>             | Polyclonal IgG                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                |