



Product Datasheet

p53 antibody (orb178524)


Description	Rabbit polyclonal antibody to p53
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human p53 (NP_000537.3).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	53KDa
Uniprot ID	P04637
Entrez	7157
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IP: 1:20-50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ma...



WB analysis of extracts of various cells...



Western blot analysis of extracts of var...