


## Product Datasheet

### SLAMF7 antibody (orb178543)

|                            |                                                          |
|----------------------------|----------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to SLAMF7                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                        |
| <b>Conjugation</b>         | Unconjugated                                             |
| <b>Tested Applications</b> | WB                                                       |
| <b>Immunogen</b>           | Recombinant fusion protein of human SLAMF7(NP_067004.3). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.         |
| <b>Concentration</b>       | 1 mg/ml                                                  |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.              |
| <b>Note</b>                | For research use only                                    |
| <b>Isotype</b>             | IgG                                                      |
| <b>Clonality</b>           | Polyclonal                                               |
| <b>MW</b>                  | 25-80KDa                                                 |
| <b>Uniprot ID</b>          | <b>Q9NQ25</b>                                            |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                          |



Western blot analysis of extracts of var...