



Product Datasheet

CD173 Antibody (PE) (orb1786025)

Description	CD173 (blood group antigen H2) is a fucosylated saccharide (Fuc-alpha-1-2-Gal-beta-1-4-GlcNAc-beta) ...
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	Human thrombocytes
Target	CD173
Preservatives	Stabilizing Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgM
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.



Flow
cytometry
surface
staining
pattern ...



Separation of
human
CD173
positive
blood...